


Termometer Telinga Inframerah / Infrared Ear Thermometer

Model: RA600



Product Reg. No.: GB63548338017

Pelindung Muncung "Plug & Secure" guna berkali / Multiple- Use PLUG & SECURE Probe Cover



Gunakan tisu alkohol untuk membersihkan pelindung muncung selepas setiap penggunaan. Gantikan pelindung muncung baru setelah 30 kali penggunaan atau bila diperlukan.
Use alcohol swab to clean the probe cover after each use. Replace with new probe cover after approximately 30 times of use or whenever necessary.

Disahkan Secara Klinikal / Clinically Validated


EN 12470-5:2003

Termometer telinga ini telah diuji secara klinikal di beberapa hospital pengajar mengikut protokol piawai peraturan EN 12470-5:2003. Ujian lengkap dengan kebolehhulangan klinikal ini dilakukan ke atas subjek yang demam dan yang mempunyai suhu badan normal.

This Ear thermometer has been clinically tested in several big teaching hospitals based on the protocol of EN12470-5:2003 regulatory standard, covering sufficient febrile and normal body temperature subjects with satisfied clinical repeatability.

- Mod Pengimbasan Berterusan 3 saat atau pengukuran 1 saat /
3 second Continuous Scanning Mode or 1 second Measurement
- Pelindung Muncung "Plug & Secure" guna berkali /
Multiple- Use PLUG & SECURE Probe Cover
- Penggera Demam /
Fever Alarm
- 9 Memori /
9 Memories
- Unit °C / °F boleh diubah /
°C / °F switchable
- Lampu Latar /
Backlight
- Padam Automatik /
Auto shut-off
- Indikasi Bateri Lemah /
Low Battery Indicator
- Diagnosis-diri untuk kerosakan /
Self-diagnosis for malfunction

HASIL DALAM 1" ATAU 3" Mod Pengimbasan RESULT IN 1" OR 3" Scanning Mode



Pelindung Muncung "Plug & Secure" guna berkali /
Multiple- Use PLUG & SECURE Probe Cover



Penggera Demam /
Fever Alarm



9 Memori /
9 Memories



Unit °C / °F boleh diubah /
°C / °F switchable



Lampu Latar /
Backlight



Padam Automatik /
Auto shut-off

Rossmax Swiss GmbH
Widnauerstrasse 1, CH-9435
Heerbrugg, Switzerland

CMSB_Epaper_RA600_ver2021