

The *Excellent Guide*

Towards Patient Care



LKL ADVANCE METALTECH SDN BHD
(Wholly-Owned Subsidiary of LKL INTERNATIONAL BERHAD)

Head Office

Wisma LKL, No. 3, Jalan BS 7/18, Kawasan Perindustrian Bukit Serdang,
43300 Seri Kembangan, Selangor Darul Ehsan, Malaysia.
Tel: +603 8948 2990 Fax: +603 8948 7904
Email: info@lklbeds.com URL: www.lklbeds.com

Penang Branch

E-32-1 Bay Avenue, Lorong Bayan Indah 2,
11900 Bayan Lepas, Pulau Pinang.
Tel: +604 640 9187/9188 Fax: +604 640 9189
Email: services_penang@lklbeds.com

Perak Branch

No 17 (Ground Floor), Laluan Tasek Timur 3,
Pusat Perdagangan Tasek Indra,
31400 Ipoh, Perak Darul Ridzuan.
Tel: +605 541 7063 Email: laiyong@lklbeds.com

Scan this QR code
for more information



Achievements & Awards

FDA Registered Company
Registration No: 3011621895



www.lklbeds.com



BA9900-WA3
HOSPITAL ICU/CCU ELECTRICAL BED

All rights reserved by LKL Advance Metaltech Sdn Bhd *copyright 2018

Due to technological and design innovations, specifications of products may be subjected to change without prior notice


BA9900-WA3

HOSPITAL ICU/CCU ELECTRICAL BED



Integrated with Weighing Scale
& Bed Exit Alarm System

Approx Dim(mm):



250kg UDL



Hi-Lo
488 - 788mm



Backrest
0° - 70°



Legrest
0° - 35°



Kneerest
0° - 20°



Trendelenburg
0° - 15°



Reverse
Trendelenburg
0° - 15°



Cardiac Chair
Position



Vascular
Position



Auto Contour



Control Buttons Embedded
Backrest Side Rail
(Inner and outer sides)




Caregiver Control Panel
with LCD Screen
Embedded Foot Board



Under Bed Light



Backrest with
X-Ray Translucent Base



Quick CPR Release Lever &
Drainage Bag Holder



Tente Twinwheel
with Central Locking



S/S IV Pole (BX305-1)

Product Features

- Zero scale function
- Bed exit alarm system and under bed light
- Nurse lock located at foot board which to lock out total or individual function of control buttons
- Quick CPR release lever for flattening the backrest located at both sides of the bed
- Auto contour function to avoid patients slipping down while performing cardiac chair position
- Hi-Lo, Backrest, Kneerest, Trendelenburg and Reverse Trendelenburg by electrical motorized
- Backrest with X-Ray Radio Translucent base complete with X-Ray cassette tray
- Four high-quality 150mm diameter German Twinwheel with central locking system (free, steer and total brakes)
- PE Head & Foot boards with locking mechanism for easy attachment, detachment and can be cleaned easily
- Split fold down PE safety side rails with one-hand operation
- Perforated Mattress support base
- Steel bed frame with scratch and chip resistant epoxy powder coating
- 6 provisions for IV pole attachment, 4 provisions for monkey pole and 2 provisions for drainage bag holders
- PU bumper at 4 corners
- IV drip pole storage

Electrical System

- Linak DC-actuator system (TUV & UL approval & EC60601-1)
- 24 VDC with built in back-up battery
- Input Voltage: 220 - 240V; 50Hz - 60Hz
- Safety Rating Control box: IPX4, IPX6 Actuators: IPX4, Nurse Control: IPX6

Accessories Included

- Pressure reducing mattress x 1 unit (BX415)
- Mattress Approx Dim (mm): 1980 (L) x 914 (W) x 125 (H)
- Telescopic IV Drip Pole with 2 hooks (Stainless steel) x 1 unit (BX305-1)

Optional Accessories




• Monkey pulley (BX910)



• Oxygen Cylinder Holder
(CWM0233)



• Portable monitor/
Defibrillator tray
(C-DEFIB-08)



• Orthopedic Frame